

LAB

EASY-CARD

FECAL SAMPLE

High definition
in gastroenterology

EASY-CARD H. pylori Ag

EASY-CARD Occult Blood

EASY-CARD Calprotectin

EASY-CARD Calprotectin Quant

SENTINEL
DIAGNOSTICS

EASY-CARD

High definition in gastroenterology

Sentinel Diagnostics Easy-Card line:

A complete profile for the definition of gastrointestinal disorders with considerable advantages for the patient and the laboratory:

- Non invasive methods: fecal sample
- Simple test execution
- Only 4 sample drops to obtain a result
- Result in 5 or 10 minutes
- Clear results interpretation
- High diagnostic reliability
- Immuno-chromatographic methods on card
- Monoclonal antibodies
- Control line for confirmation or failure

Easy-Card H. pylori Ag

Monoclonal antibody vs. H. pylori antigen.



Easy-Card Occult Blood

Immunological method, highly specific for "intact" hemoglobin in feces caused by colorectal blood loss.



Easy-Card Calprotectin

Sensitivity 50 µg/g feces.



Easy-Card Calprotectin Quant

Combination of two reaction membranes at different sensitivity: 50 and 200 µg/g feces.



Ref.	Description	Kit size
14656	Easy-Card H. pylori Ag	20 Card - 20 Collection Tube
14800	Easy-Card Occult Blood	25 Card - 25 Collection Tube
14660	Easy-Card Calprotectin	20 Card - 20 Collection Tube
14661	Easy-Card Calprotectin Quant	10 Card - 10 Collection Tube

Sentinel CH. SpA

Via Robert Koch, 2
20152 Milan - Italy
Phone: +39 02 3455141
Fax: +39 02 34551464
sentinel@sentinel.it
www.sentinel diagnostics.com



WATCHING OVER LIFE