

STATONAT[®]

LAB

Leishmania spp.

**Real-Time PCR kit
for a specific and easy
Leishmaniasis diagnosis**



SENTINEL
DIAGNOSTICS

STAT-NAT®


Leishmania spp.

LAB

For a specific and easy Leishmaniasis diagnosis.

Leishmania parasitosis is endemic in many countries and re-emerging in the Mediterranean area. A quick and accurate diagnosis allows early treatment of disease and may reduce the duration of the infection and prevent possible complications.

The **STAT-NAT® Leishmania spp.** kit, realized with STAT-NAT® technology patented by Sentinel Diagnostics, is able to detect all Leishmania species. The kit consists of a **freeze-dried master mix** that allows **room temperature transport and storage**. The STAT-NAT® Leishmania spp. kit contains all the necessary reagents for Real-Time PCR analysis.



Patent No. W02010133628 A1

Ref.	Description	Kit size
1N053	STAT-NAT® Leishmania spp.	48 test



Easy

The enzyme (Hot Start Polymerase), primers and probes are already included. One vial/one test.



Room temperature storage

STAT-NAT® is freeze-dried and guarantees long term storage at room temperature.



Performance improvement

STAT-NAT® is a ready-to-use product that minimizes analytical variables.



Whole blood

Suitable test on extracted DNA from whole blood samples.



Internal control

Amplification and extraction internal control.



Environmentally sustainable

STAT-NAT® needs neither cumbersome packaging nor temperature controlled shipment.

Sentinel CH. SpA

Via Robert Koch, 2
20152 Milan - Italy
Phone: +39 02 3455141
Fax: +39 02 34551464
sentinel@sentinel.it
www.sentinel diagnostics.com



WATCHING OVER LIFE