


**STATONAT<sup>®</sup>**

# RNA Mix

**Universal RNA Master Mix:  
an easy way to perform  
a One-Step RT-PCR**



**SENTINEL**  
DIAGNOSTICS


# STAT-NAT<sup>®</sup>

## RNA Mix

Universal RNA Master Mix: an easy way to perform a One-Step RT-PCR

The **STAT-NAT<sup>®</sup> RNA-Mix** kit, realized with STAT-NAT<sup>®</sup> technology patented by Sentinel, consists of a **freeze-dried master mix** that allows **room temperature transport and storage**.


The STAT-NAT<sup>®</sup> RNA-Mix kit is able to perform retrotranscription and nucleic acid amplification in the same vial.




Patent No.  
W02010133628 B1

The STAT-NAT<sup>®</sup> RNA-Mix kit contains all the necessary reagents for a One-Step RT-PCR reaction with the exception of primers and probes.

End Point PCR



Real-Time PCR



Ref.	Description	Kit size
1N005	STAT-NAT <sup>®</sup> RNA Mix	96 test



### Room temperature storage

STAT-NAT<sup>®</sup> is freeze-dried and guarantees long term storage at room temperature.



### Easy

STAT-NAT<sup>®</sup> contains all the reaction components. The enzymes (Hot Start Polymerase and Reverse Transcriptase) are already included.



### Performance improvement

STAT-NAT<sup>®</sup> is a ready-to-use product that minimizes analytical variables.



### Universal

STAT-NAT<sup>®</sup> technology yields very good performances in all the most diffused molecular biology techniques.



### Flexible

STAT-NAT<sup>®</sup> can improve test sensitivity even with low concentrated DNA samples.



### Environmentally sustainable

STAT-NAT<sup>®</sup> needs neither cumbersome packaging nor temperature controlled shipment.

## Sentinel CH. SpA

Via Robert Koch, 2  
20152 Milan - Italy  
Phone: +39 02 3455141  
Fax: +39 02 34551464  
sentinel@sentinel.it  
www.sentinel diagnostics.com



# SENTINEL DIAGNOSTICS



WATCHING OVER LIFE