

LAB

# EASY CARD<sup>®</sup>

FECAL SAMPLE

High definition  
in gastroenterology

**EASY CARD<sup>®</sup>** H. pylori Ag

**EASY CARD<sup>®</sup>** Occult Blood

**SENTINEL**  
DIAGNOSTICS

# EASY CARD®

## High definition in gastroenterology

### Sentinel Diagnostics EASY CARD® line:

A complete profile for the definition of gastrointestinal disorders with considerable advantages for the patient and the laboratory:

- Non invasive methods: fecal sample
- Simple test execution
- Only 4 sample drops to obtain a result
- Result in 5 or 10 minutes
- Clear results interpretation
- High diagnostic reliability
- Immuno-chromatographic methods on card
- Monoclonal antibodies
- Control line for confirmation or failure

#### EASY CARD® H. pylori Ag



Monoclonal antibody vs. H. pylori antigen.

#### EASY CARD® Occult Blood



Immunological method, highly specific for "intact" hemoglobin in feces caused by colorectal blood loss.

The test detects hemoglobin concentrations (Hb)  $\geq 0.04 \mu\text{g/mL}$  (40 ng/mL).

Ref.	Description	Kit size
14656	EASY CARD® H. pylori Ag	20 Card - 20 Collection Tube
14800	EASY CARD® Occult Blood	25 Card - 25 Collection Tube

### Sentinel CH. SpA

Via Robert Koch, 2 • 20152 Milan - Italy  
 Phone: +39 02 34 55 141 • Fax: +39 02 34 551 464  
 marketing@sentinel.it • www.sentinel diagnostics.com

