

Available
soon

FOBGold®
LINE

SENTiFIT® 800

High throughput system
for Fecal Immunochemical Testing

SENTINEL
DIAGNOSTICS

High throughput system for Fecal Immunochemical Testing


SENTiFIT® 800 Analyzer is the new high throughput fully automated system for big Fecal Immunochemical Test (FIT) routines able to process thousands samples every day. Connecting the new SENTiFIT® 800 to the SF800 Rack Loader module can improve the walk-away time and simplify the loading of a high number of samples.

FOBGold® is an immunodiagnostic (FIT) test for providing sensitive, accurate and reproducible measurements of human Hb levels in feces specimens. With its unique "2 in 1" sample collection tube, it represents the ideal choice for CRC screening programs, as well as for occult blood routine measurement. FOB Gold® assay, together with SENTiFIT® 800 Analyzer and automatic SF800 Rack Loader, is the ultimate solution to meet the needs of high volume Laboratories and Screening programs.

System specification:

| | | |
|--------------------------|------------------|--|
| SENTiFIT® 800 Analyzer | System type | Fully automated, random access and walk away instrument with piercing capability |
| | Assay type | Fecal Immunochemical Test |
| | Method | Immunturbidimetry |
| | Throughput | Up to 600 test/h |
| | Reagent capacity | Multiple reagent vials for long working sessions |
| SF800 Rack Loader Module | Sample capacity | Up to 150 samples in continuous loading |
| | Sample loading | Tubes on racks |

Specifications subject to change without notice.



| Ref. | Description |
|--------|---|
| 550800 | SENTiFIT® 800 Analyzer High throughput system for Fecal Immunochemical Testing |
| 550815 | SF800 Rack Loader Module Automatic Rack Loader for SENTiFIT® 800 Analyzer |

Sentinel CH. SpA

Via Robert Koch, 2 • 20152 Milan - Italy
Phone: +39 02 3455141 • Fax: +39 02 34551464
sentinel@sentinel.it • www.sentinel diagnostics.com

